

Alberta Aids to Daily Living

Bulletin #128

Consumer product recall

Voluntary product recall - Drive home bed assist handles recall

Alberta Aids to Daily Living (AADL)-approved manufacturer Drive DeVilbiss Healthcare (Drive), in cooperation with the U.S. Consumer Product Safety Commission and Health Canada, has informed AADL of a recent voluntary recall of certain products. The recall involves certain bed assist handles and a bed assist rail, which Drive distributed from October 2007 through June 2021.

This recall is due to serious hazard risk or death when attached to an adult's bed. Users can become trapped within the bed assist handle or rail or between the handle/rail and the side of the mattress. This poses a serious entrapment hazard and risk of death by asphyxiation.

This recall involves the four Drive models listed in the chart and pictured below. The manufacture name, Drive, and model number are printed on a label located on the product's metal tubing.

| Model | Product |
|--------------|------------------------------------|
| 15064 | Bed Assist Handle |
| 15062 | Bed Assist Rail with Folding Board |
| RTL15073 | Home Bed Assist Handle |
| RTL15063-ADJ | Home Bed Assist Handle |



15064 Bed Assist Handle



15062 Bed Assist Rail with Folding Board

December 10, 2021



RTL15063-ADJ Home Bed Assist Handle



RTL15073 Home Bed Assist Handle

Remediation plan

Each customer (dealer/distributor/retailer/re-seller) that has purchased any of the four Drive bed rails in question has received a safety recall email. The email content provides the manufacturer's official statement, remediation plan and actions customers (dealer/distributor/retailer/re-seller) need to take. The manufacturer has requested that customers destroy remaining inventory of the affected recalled products.

AADL vendors will receive similar notification of the product recall of affected products. Vendors are to immediately stop selling the recalled items. Any retailer selling, offering for sale or distributing any product covered by the recall is violating federal law.

AADL authorizers, vendors and other stakeholders will be informed of any updates via an AADL bulletin.

Action to be taken by affected clients

Immediately **stop** using the product and contact the AADL-approved vendor list, available at [AADL Mattress Overlays, Transfer Aids and Accessories](#).

More information

Jennifer Mabugat Chan, AADL Program Manager, Mobility Benefits

T: 780-422-6567

E: JenniferMabugatChan@gov.ab.ca

For more information from the manufacturer, please call (877) 467-3099 from 8:00 a.m. to 7:00 p.m. ET, Monday through Friday.